

Microplate Enzyme Assay Using Fluorescence

Original Reference: [Include a citable reference even if you got a protocol from another lab or this website]

Protocol Prepared By: Meg Steinweg (steinweg@nrel.colostate.edu) & Shawna McMahon (shawna.mcmahon@gmail.com)

Date Created: July 18, 2007

Last Updated: May 13, 2009

This assay is applicable to any enzyme for which there is a fluorescently-labelled (MUB or MUC) substrate:

- β -D-cellubiosidase (CB)
- α -Glucosidase (AG)
- β -Glucosidase (BG)
- L-Leucine-7-amidomethylcoumarin (LAP)
- N-acetyl- β -glucoasminidase (NAG)
- Phosphatase (PHOS)
- β -Xylosidase (XYL)

Equipment

- 15 mL centrifuge tubes
- Pipetting reservoirs (one for each standard concentration (6) x two standards = 12 + one for each substrate = 18 for all 6 substrates)
- Centrifuge equipped with holders for deep-well plates
- Deep-well plates with lids (8 plates for 12 samples—4 each for samples & standards at 4 temps.)
- Waring Blender with small stainless steel blending cup
- Small glass dishes
- Magnetic stir bars and stir plate
- 8-channel repeat pipettor (volume distribution at least 200uL-800uL)

Reagents & Solutions

- 4-MUB-B-D- xylopyranoside: Sigma M7008 (100mg)
- 4-MUB- β -D-glucopyranoside: Sigma M3633 (250mg)
- 4-MUB- α -D-glucopyranoside: Sigma M9766 (100mg)
- 4-MUB-N-acetyl- β -D-glucosaminide: Sigma M2133 (100mg)
- 4-MUB- β -D-cellobioside: Sigma M6018 (100mg)
- 4-MUB phosphate: Sigma A8625 (25g)
- L-leucine-7-amido-4-methylcoumarin hydrochloride: Sigma L2145 (25 mg)
- 7-amino-4-methylcoumarin (MUC): Sigma A9891 (250mg)
- 4-methylumbelliferone (MUB): Sigma M1381
- 50 mM Buffer
 - * Buffer pH should match the soil pH as closely as possible
 - for acidic soils → sodium acetate buffer
 - for basic soils → Tris buffer
 - * <http://www.bioinformatics.org/JaMBW/5/4/index.html> is a good buffer calculator

Fluorescent Substrate Preparation (200 μ M)

Dissolve the fluorescent substrates in deionized water according to the following list:

- β -Xylosidase4-MUB- β -D- xylopyranoside: 6.17 mg/100 mL DI H₂O
- BG: 4-MUB- β -D-glucopyranoside: 6.77 mg/100 mL DI H₂O
- AG: 4-MUB- α -D-glucopyranoside: 6.77 mg/100 mL DI H₂O
- NAG: 4-MUB-N-acetyl- β -D-glucosaminide: 7.59 mg/100 mL DI H₂O
- CB: 4-MUB- β -D-cellobioside: 10 mg/100 mL DI H₂O
- Phosphatase: 4-MUB phosphate: 5.12 mg/100 mL DI H₂O
- LAP: L-leucine-7-amido-4-methylcoumarin hydrochloride: 6.5 mg/100 mL DI H₂O

Stir and use moderate heat to speed dissolution of substrates.

Aliquots: Put 11 mL aliquots into 15 mL centrifuge tubes.

Storage: Can be stored frozen (-20°C for up to two months). Thaw in warm water. Do not thaw and refreeze

Stock Standard Solution Preparation (1 mM)

Dissolve 17.6 mg of MUB in 100 mL DI H₂O (use a volumetric flask).

Dissolve 17.5 mg of MUC in 100 mL DI H₂O (use a volumetric flask).

Aliquots: Put 5 mL aliquots into 15 mL centrifuge tubes.

Storage: Can be stored frozen (-20°C for up to two months). Do not thaw and refreeze.

Standard Curve

- Dilute standard stock solution (1 mM=1000 μ M) to create a standard curve. Start with a 1:10 dilution of the stock solution to create a 100 μ M solution. Make all further dilutions from 100 μ M. Aim for 30 mL of each standard for a typical assay run with 12 samples.
 - 100 μ M = 3 mL 1000 μ M stock into 27 mL DI H₂O (1:10 dilution)
 - 50 μ M = 15 mL 100 μ M into 15 mL DI H₂O (1:2 dilution)
 - 25 μ M = 15 mL 50 μ M into 15 mL DI H₂O (1:2 dilution)
 - 10 μ M = 12 mL 25 μ M into 18 mL DI H₂O (1:10 dilution)
 - 5 μ M = 15 mL 10 μ M into 15 mL DI H₂O (1:2 dilution)
 - 2.5 μ M = 15 mL 5 μ M into 15 mL DI H₂O (1:2 dilution)
- To make different amounts of standard, use the following formula:

$$(\text{High Concentration})(x) = (\text{Low Concentration})(\text{Desired Volume})$$

x = volume of the higher concentration to add to water

add to volume of water = (desired volume – volume of higher concentration)

e.g. To make 80 mL of 25 μ M from 50 μ M:

$$(50 \mu\text{M})(x) = (25 \mu\text{M})(80 \text{ mL})$$

$$x = 40 \text{ mL of } 50 \mu\text{M in } 40 \text{ mL H}_2\text{O}$$

- Because of background autofluorescence, a standard curve must be prepared for each sample. See plate layout for clarification.

Buffer Preparation

- The buffer you select will depend on the pH of the soil. Below are the instructions for two common buffers: sodium acetate (pKa = 4.76) and Tris (pKa = 8.06).
 - Note that phosphate buffer has a pKa = 7.2 which could work for more neutral/slightly basic soils. However, high phosphate concentrations may interfere with enzyme activity so use with caution.

50 mM Sodium Acetate Buffer

1. Dissolve 30.8 g sodium acetate in 2L DI H₂O
2. Adjust the pH to match soil pH using glacial acetic acid or concentrated sodium hydroxide (~10M NaOH).
3. Bring the volume up to 4L with DI H₂O.

Storage: Buffer should be made fresh, but can be refrigerated (4°C) for up to a week.

50 mM Tris Buffer

- * notes on this are in my notebook – will add later!

Logistics

- Standards and substrates should be prepared well ahead of time.
- On the morning of the assay:
 - Thaw substrates and required standards (MUC only needs to be thawed if you are doing the LAP assay).
 - Prepare buffer.
 - Prepare standard curve (see below for dilutions)
- This method is best done with two people (one preparing soil slurries and the other preparing plates). Two people can easily prepare assays for up to **12 soil samples** in a morning (at most 2 hrs).
- Expert practitioners can start two sequential sets of assays in a morning. If you start assays in the afternoon, you will likely have to come in late to read the 15°C plates. Be sure to note incubation times for temperature of interest:

Temperature	Time
4°C	~ 23 hours
15°C	6 hours
25°C	3 hours
35°C	1.5 hours

- Preparing samples to be read on the microplate reader takes one person approximately 10 min. per plate. Total analysis time is typically 1.5 hours over the course of a day for a full set of samples.

Tips & Tricks

- You may need to alter the plate layout for your experiment. However, we've provided the plate layouts we use for doing all the enzymes at once.
- Because of differences in background autofluorescence, you must **prepare a standard curve for each sample** separately, hence the large standard plates.

Plate Layouts:

Enzyme Activity Samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	S1+ BG	S2+ BG	S3+ BG	S4+ BG	S5+ BG	S6+ BG	S7+ BG	S8+ BG	S9+ BG	S10+ BG	S11+ BG	S12+ BG
B	S1+CB	S2+CB	S3+CB	S4+CB	S5+CB	S6+CB	S7+CB	S8+CB	S9+CB	S10+CB	S11+CB	S12+CB
C	S1+ NAG	S2+ NAG	S3+NAG	S4+ NAG	S5+ NAG	S6+ NAG	S7+ NAG	S8+ NAG	S9+ NAG	S10+ NAG	S11+ NAG	S12+ NAG
D	S1+PHOS	S2+PHOS	S3+PHOS	S4+PHOS	S5+PHOS	S6+PHOS	S7+PHOS	S8+PHOS	S9+PHOS	S10+PHOS	S11+PHOS	S12+PHOS
E	S1+ XYL	S2+XYL	S3+ XYL	S4+XYL	S5+XYL	S6+XYL	S7+XYL	S8+XYL	S9+XYL	S10+XYL	S11+XYL	S12+XYL
F												
G												
H	S1 +LAP	S2 +LAP	S3 +LAP	S4+LAP	S5+LAP	S6+LAP	S7+LAP	S8+LAP	S9+LAP	S10+LAP	S11+LAP	S12+LAP

4 and 25C Standards Same for MUB and MUC

	1	2	3	4	5	6	7	8	9	10	11	12
A	S1+ std0	S2+ std0	S3+ std0	S4+ std0	S5+ std0	S6+ std0	S7+ std0	S8+ std0	S9+ std0	S10+ std0	S11+ std0	S12+ std0
B	S1+std2.5	S2+std2.5	S3+std2.5	S4+std2.5	S5+std2.5	S6+std2.5	S7+std2.5	S8+std2.5	S9+std2.5	S10+std2.5	S11+std2.5	S12+std2.5
C	S1+ std5	S2+ std5	S3+ std5	S4+ std5	S5+ std5	S6+ std5	S7+ std5	S8+ std5	S9+ std5	S10+ std5	S11+ std5	S12+ std5
D	S1+std10	S2+std10	S3+std10	S4+std10	S5+std10	S6+std10	S7+std10	S8+std10	S9+std10	S10+std10	S11+std10	S12+std10
E	S1+ std25	S2+ std25	S3+ std25	S4+ std25	S5+ std25	S6+ std25	S7+ std25	S8+ std25	S9+ std25	S10+ std25	S11+ std25	S12+ std25
F	S1+ std50	S2+ std50	S3+ std50	S4+ std50	S5+ std50	S6+ std50	S7+ std50	S8+ std50	S9+ std50	S10+ std50	S11+ std50	S12+ std50
G	S1+ std100	S2+ std100	S3+ std100	S4+ std100	S5+ std100	S6+ std100	S7+ std100	S8+ std100	S9+ std100	S10+ std100	S11+ std100	S12+ std100
H												

- The spaces in the enzyme activity sample plate at rows F and G are there in case other MUB based enzymes need to be measured at some point.
- S1...S12 = samples 1-12
- std0...std100 = standard concentrations

Procedure

I. Assay Set-Up

1. Prepare standard plates by pipetting 200 μL of appropriate standards into correct wells of MUB or MUC standard plates (oriented in *rows*).
 - * When you add soil in Step 5, record the time and read these plates after the same amount of elapsed time as the enzyme assay plates at the corresponding temperature.
 - Prepare two plates for each standard (MUB & MUC) for a total of four plates. One set is incubated at 4°C and the other at 25°C.
 - The 25°C standard curve is used for all but the 4°C assay (15°C, 25°C & 35°C).
2. Prepare soil slurries by weighing 2.75 grams of field moist soil into blender.
 - The soil amount can be altered, just be sure to adjust the buffer amount proportionally.
 - Be sure to measure the gravimetric water content so that dry weight (d.w.) can be calculated—enzyme activity is typically presented on a g^{-1} soil d.w. basis.
3. Add 91 mL of 50 mM buffer.
4. Blend contents on high for 1 min.
5. Pour contents into a glass bowl, add stir bar and place on stir plate. Mix on low.
 - Bowl must be wide enough to accommodate the 8-channel pipette.
 - Stirring gently helps keep soil suspended and minimize variation due to different amounts of soil present in each well.
6. Pipette 800 μL of soil slurry into wells for enzyme activity measurement and standard curves (see plate layouts). Each sample should be in a separate *column*.
7. Rinse blender with DI H_2O or buffer between samples. Repeat Steps 2-7 for each sample.
8. Pipette 200 μL of appropriate 200 μM substrate into correct assays wells (oriented in *rows*).
 - * Record the time of substrate addition and try to start all plates/substrates as quickly as possible.

Note: Steps 1 & 8 are mostly easily performed with only six (6) pipette tips loaded on the repeat pipette (two additions across the plate for a total of 12).

II. Incubation & Analysis

1. Seal the plates with plate mats and mix by hand (invert about 10 times).
2. Place plates in appropriate incubators for required incubation period.

Temperature	Time
4°C	~ 23 hours
15°C	6 hours
25°C	3 hours
35°C	1.5 hours
3. When incubation is complete, centrifuge plates for 3 min. at 1500 rpm (~350 $\times g$).
4. Transfer 250 μL from each well into corresponding well in a flat-bottomed black 96-well plate.
5. Add 5 μL of 0.5N NaOH to each sample.
 - * Want final pH in wells to be ~10.
 - Adjust normality as necessary for pH of soil.

Note: this step may be omitted without significant loss of sensitivity in high pH soils.

6. Following manufacturer's instructions for your plate reader, measure fluorescence using the following parameters:

Excitation Wavelength = 365

Emission Wavelength = 450

Additional instructions for the Tecan Infinite M200 microplate reader:

1. Select 10 reads but no multiple reads per well.
2. Open Excel and make sure the correct destination file is open.
 - The Tecan will put the data in a new worksheet (tab) in whatever spreadsheet is open. Label each tab or put descriptive information in each worksheet as it is generated. Information needed, sample #s, MUB or MUC sample, MUB or MUC standard, gain of the microplate reader. If you are using numbers that don't mean anything, then do a summary sheet on each workbook that explains the numbers. e.g. Sample 1 is plot 25 0-5cm, Sample 2 is plot 25 5-15cm, etc...
2. Run the standard plate at optimal gain. The highest measured value should be for the 100 μM standard. Note the optimal gain determined by the machine. If you are doing the 35°C temp then you won't have the standard information yet, so be careful, you may need to drop the gain more than 5 units. I would suggest dropping it 5 units, run the sample and then drop it again 5 units and run the sample again. That way when you run the 25°C standards 1.5hr later you have another option if the 100 μM standard is still too high for the first 5unit gain drop (another trick of the trade).
3. Set the gain manually to 5 units below the optimal gain from the standard plate for all other plates with the same soil samples.
 - * The gain dramatically changes the measured values. However, it may be necessary to drop the gain even further for some samples if the value exceeds detection limits (indicated by a red OVER in the spreadsheet output). It is **CRITICAL** that the standard plates are read at the same gain as the samples else the sample values are **MEANINGLESS**.

Data Analysis & Calculations

- Plot standard curves (fluorescence on Y, standard conc. on X) in Excel for the MUB and MUC standards for each soil sample. Calculate the slope, intercept & r^2 . Accept the standard curve if r^2 is >0.98 .
 - We often do two standard curves for each set of standards to get better resolution for high activity and low activity enzymes, but this isn't strictly necessary as the standard curves are quite linear.
 - Low Activity: plot 0-25 μM standards
 - High Activity: plot 10-100 μM standards
- Use $y=mx+b$ to solve for unknown concentrations (x), where:
 - y = sample fluorescence
 - x = enzyme concentration
 - m = slope from standard curve
 - b = intercept from standard curve

sample enzyme concentration = (sample fluorescence – std. curve intercept) \div std. curve slope

- Calculating the enzyme activity:

Enzyme Concentration/Length of Incubation (h)/soil dry weight (g) = $\mu\text{mol activity/g soil d.w./h}$

Notes

Safety

- This is a relatively benign analytical technique. Wear gloves if you are sensitive and consider safety glasses.